LABORATORY REQUISITION FORM FOR HEPATITIS C RESISTANCE TEST

A) PATIENT, REQUESTING CLINICIAN, AND REFERRING LABORATORY	
Patient's Last	BC Ordering Clinician
Name: Patient's First Name:	Clinician in other provinces: Submitter's name should be the Public Health or Regional Laboratory Physician who will be specified by that lab in the space on the left
Date of Birth: DD MMM YYYY	Clinician's Name & MSP No.: Street: City, Province, Postal Code: Phone: Fax: Sample Comment Cc Dr/NP Name
Province: please circle PLNA PLSA MB SK ON NB NL NS PE YT	& MSP No.:
B) HCV RESISTANCE TEST	
This test is recommended for patients who are living with Her particular, those who have failed a previous HCV therapy conconsidered as being for <i>Research Use Only</i> . Patient's HCV Genotype (if available):	
Plasma sample was already referred to BCCDC PHL Virology for other HCV tests? Check if yes.	DD MMM YYYY
This space for BC CfE use only: Sample info provided to BCCDC PHL Viro on date:	
Plasma sample is being shipped to BC-CfE Lab for this test? Check if yes. This space beside is for BC CfE use only: Sample ID label	Write any other specimen collection dates, in the same format as above, if additional samples from the same patient are desired for testing
Request Date:	
	Name of requestor if other than the ordering clinician or referring physician:
C) SPECIMEN COLLECTION, STORAGE, AND TRANSPORT INSTRUCTIONS	

To BC Labs: Our lab has access to the patient's plasma from HCV viral load testing performed at the BCCDC PHL, so there is no need to send us another sample. *Providing the specimen collection date above* will be sufficient to submit this test request by faxing it to us.

Specimen Collection: A minimum of 1.2 mL of EDTA **PLASMA** is required. Collect one appropriate volume EDTA vial. Centrifuge EDTA vial for 15 mins at 800-1600 rcf. Transfer >1.2 mL of plasma to one 2 mL screwcap cryovial.

Storage: Plasma can be stored for up to 48 hrs at 2° to 8° C. Freeze at -15° to -80° C during delays to our lab of >48 hrs after collection.

Transport: For local Metro Vancouver shipments of refrigerated plasma vials, ship on/between a suitable number of **well-chilled or frozen ice packs.**

From labs outside of Metro Vancouver and outside of BC: Ship frozen EDTA plasma on dry ice. Ship by overnight carrier on Monday to Wednesday only, allowing 2 consecutive regular business days for delivery in case of any weekday transport delays.

Notify our lab of your shipment (fax, phone, or email - but do not include patient identification by email), include your **carrier company and its tracking number**, note how many samples and which test types if for more than HCV. Follow IATA and TDG Regulations (UN3373, Category B).

Ship to: BC Centre for Excellence in HIV/AIDS, Research Laboratory Room 604, Attn: Dr. Chanson Brumme, St. Paul's Hospital, 1081 Burrard Street, Vancouver, BC V6Z 1Y6 | Lab Business Tel: 604.806.86775 | lab clerks@bccfe.ca

BC-CfE Research Laboratory test results, like all laboratory results, are kept strictly confidential.

Further information and order forms can be found at http://bccfe.ca/laboratory-program/laboratory-test-order-forms